



One Step
ahead **NOVO**

Novo Klinik-Service GmbH

Product Catalogue

At Your Service

Tracheostomy / Laryngectomy
fixations and accessories

Ventilation
fixations and accessories

Securement Devices and Prophylaxis

Customised Products

About Us

You About Us

References

Contact

Novo Klinik-Service ... One Step ahead

Since 1987 **NOVO** has specialised in the development and production of latex-free fixations, cushions and medical devices, all for the ventilated patient.

Our products are widely spread in clinics as well as in the care, home care and rescue sector.

Besides safety our aid devices offer the highest comfort for patients and medical staff. Predominantly breathable and extremely soft materials protect the skin from irritations and decubitus ulcer. The convenient handling offers quick and safe application.

Due to extensive experience and clinical trials in collaboration with national and international institutions, we have pioneered innovative fixations. We incorporate your ideas and recommendations gained in the practical application in order to develop targeted and novel products.

If there is a problem, we will always strive to solve it.
Just contact us!

We are here for you, world wide

Austria, Bahrain, Belgium,
Canada, Costa Rica,
Czech Republic, Denmark,
El Salvador, Estonia, Finland,
France, Great Britain, Greece,
Hungary, Iceland, Ireland,
Italy, Jordan, Latvia, Macedonia,
Montenegro, Poland, Portugal,
Saudi Arabia, Serbia, Slovakia,
Slovenia, Spain, Sweden,
Switzerland, The Netherlands,
United Arab Emirates

Our global attendance on trade shows and our training see to the application of products on the patient at their best.




Personal assistance and direct customer supply go without saying for us.



! the Original !

BePa-Clip®


One-Piece Tracheostomy Tube Holder




BePa-Clip VARIO 25/K



BePa-Clip VARIO 36/H,
Novo-Pad 2



BePa-Clip VARIO 25/H,
Novo-Pad



BePa-Clip VARIO 30/H

BePa-Clip® is an one-piece fixation, which safely secures tracheostomy tubes on adults and adolescents. The slightly absorbing, breathable padding offers the best comfort with an optimum protection against decubitus ulcer and highest dimension stability.

BePa-Clip® is safe as well as quick and easy to use. It is available in different sizes and optionally with a hook or Velcro® fastener.

We can also customise the product according to your needs.


	Fastener	Sales Unit	REF	CE	RoHS	REACH	LF	LF
BePa-Clip® VARIO 18/H	hook	50 pcs	200 301					
BePa-Clip® VARIO 18/K	Velcro®	50 pcs	200 302					
BePa-Clip® VARIO 25/H	hook	60 pcs	200 501					
BePa-Clip® VARIO 25/K	Velcro®	60 pcs	200 502					
BePa-Clip® VARIO 30/H	hook	60 pcs	200 601					
BePa-Clip® VARIO 30/K	Velcro®	60 pcs	200 602					
BePa-Clip® VARIO 36/H	hook	60 pcs	200 701					
BePa-Clip® VARIO 36/K	Velcro®	60 pcs	200 702					

! the Original !




BePa-KidsClip®

One-Piece Tracheostomy Tube Holder for premature infants, babies and children



BePa-KidsClip
VARIO 13/K



BePa-KidsClip
VARIO 13/H

BePa-KidsClip®, the tracheostomy tube holder for premature infants, babies and children is particularly soft and protects the tender skin from decubitus ulcer. Offering the utmost comfort the slightly absorbing and breathable padding ensures an excellent dimension stability.

BePa-KidsClip® is safe as well as quick and easy to use. It is available in different sizes and optionally with a hook or Velcro® fastener.


We can also customise the product according to your needs. Just contact us!

	Fastener	Sales Unit	REF	CE	RoHS	REACH	LF	LF
BePa-KidsClip® VARIO 11/H	hook	60/100 pcs	200 201					
BePa-KidsClip® VARIO 11/K *	Velcro®	60/100 pcs	200 202					
BePa-KidsClip® VARIO 13/H	hook	60/100 pcs	200 211					
BePa-KidsClip® VARIO 13/K	Velcro®	60/100 pcs	200 212					
BePa-KidsClip® VARIO 15/H	hook	60/100 pcs	200 221					
BePa-KidsClip® VARIO 15/K	Velcro®	60/100 pcs	200 222					


* with extremely narrow elastic and plastic slider

Duo-Clip


Two-Part Tracheostomy Tube Holder



Duo-Clip K



Duo-Clip H



Duo-Clip H, Novo-Pad

Duo-Clip is a two-part tracheostomy tube holder for adults, adolescents and children.

This fixation represents the synergy of comfort and flexibility. The soft cotton-covered padding offers the ideal decubitus ulcer protection and is quick and easy to apply. The Velcro® fastener is securely attached to the padding by the technique of ultrasonic welding. Skin irritations caused by joining seams are consequently precluded.

The fixation is infinitely variable and can hence be adjusted to any neck girth.

Duo-Clip is optionally available with a hook or Velcro® fastener.

We can also customise the product according to your needs.


Just contact us!

	Fastener	Sales Unit	REF
Duo-Clip H	hook	60 pcs	200 801
Duo-Clip K	Velcro®	60 pcs	200 802




Duo-KidsClip

Two-Part Tracheostomy Tube Holder for babies, toddlers and children



Duo-KidsClip



Duo-KidsClip

Duo-KidsClip, the two-part tracheostomy tube holder for babies, toddlers and children offers a perfect decubitus ulcer protection for sensitive skin by its softness of the cotton covered padding.

The fixation's comfort as well as the infinite adjustment facilitate a quick and easy application to the user and the little patient.

Duo-KidsClip is an absolutely seamless holder. Due to the innovative technique of ultrasonic welding the Velcro® fastener is securely attached to the padding. Skin irritations caused by seams or overhanging parts are, therefore, impossible.


The fixation is available with Velcro® fastener.

	Fastener	Sales Unit	REF
Duo-KidsClip	Velcro®	60 pcs	200 805



Trachea-Fix

One-Piece Tracheostomy Tube Holder



Trachea-Fix, Novo-Pad ultrasonic welding

Trachea-Fix is a one-piece fixation, which safely secures tracheostomy tubes on adults and adolescents.

On the one hand this holder is preferentially applied at just operated patients, on the other hand it is especially made for long-term carriers of tracheal cannulae as it is soft, discreet and dimensionally stable at the same time.

Due to its Velcro® fastener this fixation is infinitely variable and therefore appropriate for each neck girth. The fastener is firmly attached to the padding by ultrasonic welding so that skin irritations caused by joining seams are consequently precluded.

Trachea-Fix is available with Velcro® fastener in one universal size and exceedingly quick and easy to apply.

	Fastener	Sales Unit	REF
Trachea-Fix	Velcro®	60 pcs	200 442



Sim-Plex

One-Piece Fixation for Tracheal Cannulae, Mouthguards and Oxygen Masks



Sim-Plex, Novo-Pad hook fastener

Sim-Plex is an extremely discreet one-piece fixation made from woven, soft elastic equipped with a hook fastener.

It is used to fix tracheal cannulas for adults and adolescents as well as to fix oxygen masks. In endoscopy the holder serves as fixation for the mouthguard.

Thanks to the plastic slider **Sim-Plex** is infinitely adjustable and therefore available in one universal size.

	Fastener	Sales Unit	REF
Sim-Plex	hook	60 pcs	200 445





Novo-Pad


Tracheostomy Foam Dressing




Novo-Pad, BePa-Clip VARIO 30/K



Novo-Pad, Duo-Clip H



Novo-Pad, BePa-Clip VARIO 25/H



Novo-Pad 2, BePa-Clip VARIO 36/H

Novo-Pad is a soft, skin-friendly and absorbent tracheostomy dressing for adults and adolescents.

It is made from a hydrophilic polyurethane foam that prevents the skin from pressure marks as a cushioning. Additionally, it provides an ideal removal of wound secretions with no extensive seepage through the cover.

The skin-coloured cover is semipermeable so that the absorbed moisture is gradually wicked away from the surface. Bacteria and germs cannot penetrate from the outside.

The sterile dressing **Novo-Pad** ensures an optimal wound milieu. As it does not adhere to the skin, it can be removed without trauma. The wound healing will be supported.

It is available in two sizes.

	Length	Width	Sales Unit	REF	CE 0535	UNO	STERILE	DIN EN 980
Novo-Pad	8,2 cm	9,5 cm	25 pcs	70 - 10				
Novo-Pad 2	8 cm	14 cm	25 pcs	70 - 30				

Novo-KidsPad

Tracheostomy Foam Dressing for babies, toddlers and children



Novo-KidsPad, BePa-KidsClip VARIO 13/K



Novo-KidsPad, BePa-KidsClip VARIO 13/K

Novo-KidsPad, the skin-friendly and absorbent tracheostomy dressing for children, is especially soft and offers an optimal cushioning for the tender skin.

The white layer made from a hydrophilic polyurethane foam reliably removes the wound secretion into the inside of the dressing. The skin-coloured cover is semipermeable. The absorbed moisture is gradually wicked away from the surface. Bacteria and germs cannot penetrate through the layer from the outside.

Novo-KidsPad prevents maceration of the wound. It does not adhere to the skin.

The zigzag slit facilitates application and prevents the tracheostomy dressing slipping even with mobile children.

	Length	Width	Sales Unit	REF	CE 0535	UNO	STERILE	DIN EN 980
Novo-KidsPad	6,5 cm	6,5 cm	25 pcs	70 - 20				

non-adhesive
TuBo-Grip
 Two-Part Endotracheal Tube Holder




TuBo-Grip complete extra long



TuBo-Grip complete



TuBo-Grip complete



TuBo-Grip complete

TuBo-Grip is a two-part endotracheal tube holder for adults, adolescents and pupils.

The middle part has an elastic strap with a silicone layer, which clings to the tube and offers stability even with intensely salivating patients. The neck padding is wide, softly cushioned and therefore specially convenient for long-term intubation. It is gladly used with sensitive skin and obese patients.

Due to the Velcro® fastener the fixation can be opened and readjusted at any time. Hence it is especially suitable for the field of endoscopy, too.

TuBo-Grip is optionally available in two different neck padding sizes.

	Length	Sales Unit	REF	CE	i	no	no	no	no
TuBo-Grip complete		60 pcs	20 - 50						
TuBo-Grip complete extra long		60 pcs	20 - 51						
TuBo-Grip neck padding	32 cm	60 pcs	20 - 52						
TuBo-Grip neck padding extra long	47 cm	60 pcs	20 - 53						
TuBo-Grip middle part	33 cm	60 pcs	20 - 55						

non-adhesive

Soft-Fix

Multi-Fixation



Soft-Fix normal,
Soft-Fix Velcro®



Soft-Fix small,
Soft-Fix Velcro®



Soft-Fix large,
Soft-Fix Velcro®



Soft-Fix small,
Soft-Fix Velcro®

Soft-Fix is a non-adhesive, multi-functional fixation that can be applied for patients of any age. It is used to protect and stuff infusion tubes, probes, drainages, catheters, cables and tube systems for life-maintaining technical apparatus.

It is slightly absorbing and relockable so that it is used for the fixation of wound dressings and as a decubitus prophylaxis.


Soft-Fix helps to reduce or avoid tensile loading, pressure points and decubitus ulcer.

Soft-Fix is available on rolls of three different widths or as a mixed unit each packed in a handy dispenser box.


	Width	Length	Sales Unit	REF	CE	EN	ISO	ISO	ISO	ISO
Soft-Fix small	2 cm	5 m	12 pcs	30 - 01						
Soft-Fix normal	4 cm	5 m	6 pcs	30 - 02						
Soft-Fix large	6 cm	5 m	4 pcs	30 - 03						
Soft-Fix Velcro® roll	2 cm	15 m	1 pcs	30 - 04						
Soft-Fix Mix	5 x small, 2 x normal									
	1 x large, 1 x Velcro® roll		9 pcs	30 - 05						

Main-Lock Set


Adhesive Tape for Catheters and Drainages



Main-Lock Set 14



Main-Lock Set 14



Main-Lock Set 14



Main-Lock Set 7

Main-Lock Set is a skin-friendly tape with a hook and loop connection by means of which long-term catheter and drainage tubes can be securely fixed. This leads to releasing the traction on the puncture and therefore reducing catheter-associated infections.

Main-Lock Set may be used wherever non-adhesive fixations are not desired or applied due to physiological reasons. The main field of application is the security of central venous catheters at the aditus of the subclavian vein, jugular vein as well as the femoral vein and artery. Furthermore, the adhesive tape is applied for thoracic drainages, suprapubic urinary catheters and PEG probes, in infusion therapy as well as in dialysis.

Main-Lock Set is available in two sizes.

	Width	Length	Sales Unit	REF	CE	EN	ISO	ISO	ISO	ISO
Main-Lock Set 7	4 cm	7 cm	60 pcs	90 - 30 - 10						
Main-Lock Set 14	4 cm	14 cm	40 pcs	90 - 30 - 20						

Splint-Fix

Multifunctional Splint



Splint-Fix 9



Splint-Fix 16



Splint-Fix 22



Splint-Fix 22

Splint-Fix is a softly padded one-piece splint for hands or feet serving as an immobilization and support for the wrist respectively ankle in case of arterial accesses and peripheral infusions.

The application of **Splint-Fix** ensures a continuous supply of medicine. Furthermore, it prevents the formation of a hematoma as well as unintended, paravenous infiltrates.

This splint can be individually formed according to physiological requirements and is permeable to X-rays. It can also remain at the patient during computer tomography.


It is available in five sizes.

	Width	Length	Sales Unit	REF	CE	i	RoHS	REACH	Flame Retardant	Medical Device
Splint-Fix 9	3-4 cm	9 cm	20 pcs	60 - 21						
Splint-Fix 13	4 cm	13 cm	20 pcs	60 - 22						
Splint-Fix 16	5 cm	16 cm	20 pcs	60 - 23						
Splint-Fix 22	9 cm	22 cm	14 pcs	60 - 11						
Splint-Fix 33	9 cm	33 cm	10 pcs	60 - 13						


non-adhesive 

PulsOxy-Fix

One-Piece Fixation for Pulse Oximeter Sensors



PulsOxy-Fix 17



PulsOxy-Fix 12



PulsOxy-Fix 12



PulsOxy-Fix 10

PulsOxy-Fix is a softly padded fixation for pulse oximeter sensors. It is fixed on the skin non-adhesively and can be applied at the foot, finger or other extremities.

Due to the Velcro® fastener the sleeve can be opened and closed at any time. Hence, relocating the sensor is easily and fastly done without changing the fixation.

PulsOxy-Fix is available in three sizes.

	Width	Length	Sales Unit	REF	CE	i	RoHS	REACH	Flame Retardant	Medical Device
PulsOxy-Fix 10	1,8 cm	10 cm	120 pcs	30 - 14						
PulsOxy-Fix 12	3 cm	12 cm	80 pcs	30 - 10						
PulsOxy-Fix 17	4 cm	17 cm	60 pcs	30 - 12						

Clip-Fix


One-Piece Tube Holder



Clip-Fix 1



Clip-Fix 1



Clip-Fix 3, Uni-Fix



Clip-Fix 3 DIN

Clip-Fix is a Velcro® strap with a metal holder, which serves for the sorting of tube and cable systems at the bed sheet, patient's clothing or at the DIN rail.

For emergency rescue services **Clip-Fix** is the optimal facility to sort cables, e.g. EGG cable, to fix catheter or infusion bags in the mobile ICU vehicle, but also directly at the scene and when transporting patients.


The fixation is a reusable product, it is washable at 60 °C and can be disinfected by standard disinfectant or vapour sterilisation. The metal clip consists of an acid-resistant steel.

Clip-Fix is available in three different versions and in different lengths.


	Width	Length	Sales Unit	REF	CE	IP	RoHS	REACH	EN60601	EN60601-1
Clip-Fix 1	2,5 cm	11 cm	50 pcs	30 - 30 - 01						
Clip-Fix 2	2,5 cm	17 cm	50 pcs	30 - 30 - 02						
Clip-Fix 3	2,5 cm	23,5 cm	40 pcs	30 - 30 - 03						
Clip-Fix 1 Duo	2,5 cm	12,5 cm	50 pcs	30 - 30 - 04						
Clip-Fix 2 DIN	2,5 cm	18,5 cm	50 pcs	30 - 30 - 02D						
Clip-Fix 3 DIN	2,5 cm	25 cm	40 pcs	30 - 30 - 03D						

Hand-Fix

Hand and Foot Fixation



Hand-Fix 25



Hand-Fix 30



Hand-Fix 12



Hand-Fix 17

Hand-Fix is a soft and stuffed fixation for hands and feet. It is used for self-protection of patients by preventing unintended manipulation of inlets, bandages or life-sustaining technical devices.

The used materials are skin-friendly and soft; hence they do not incise the skin.


This fixation is outstandingly appropriate to tender or sensitive skin. The sleeve and the strap ends can be cut to size.

Hand-Fix is available in five sizes.


	Width	Length	Sales Unit	REF	CE	IP	RoHS	REACH	EN60601	EN60601-1
Hand-Fix 12	3 cm	12 cm	30 pcs	60 - 01						
Hand-Fix 17	4 cm	17 cm	30 pcs	60 - 08						
Hand-Fix Mix	14 x HF/12, 14 x HF/17		28 pcs	60 - 50						
Hand-Fix 22	7 cm	22 cm	16 pcs	60 - 61						
Hand-Fix 25	7 cm	25 cm	16 pcs	60 - 63						
Hand-Fix 30	7 cm	30 cm	16 pcs	60 - 65						

Oxy-Care / Ute-Fix


Decubitus-Prophylaxis



Oxy-Care



Ute-Fix 2



Ute-Fix 3

Oxy-Care / Ute-Fix are comfortable, atraumatic paddings to protect sensitive skin. They can basically be attached to any tube system.

Oxy-Care protects particularly tender skin at the ear when applying a nasal cannula.

It is available in one size.

Ute-Fix is especially applied for probes and drainages in order to avoid respectively reduce the risk of secondary harms like ulcerisation, maceration and the potential risk of skin infection.

Ute-Fix is available in three sizes.

	Ø	Length	Sales Unit	REF	CE	i	RoHS	REACH	Star	Umbrella
Oxy-Care	1 cm	7 cm	120 pcs	40 - 10						
Ute-Fix 1	1,5 cm	2 cm	60 pcs	40 - 06						
Ute-Fix 2	2 cm	2 cm	60 pcs	40 - 05						
Ute-Fix 3	2,5 cm	2 cm	60 pcs	40 - 01						

One Step
ahead
??!

Your product?
Just contact us!

Individual solutions according to your needs

The ideal patient care also needs solutions of specific problems. This is why we develop and manufacture products customised in size and type.

Are you in need of a customised fixation? Just contact us!

BePa-KidsClip

special sizes:
e.g. 7 cm and 9 cm



TuBo-SA

customised fixation for
endotracheal tubes



BePa-Clip Standard

special size:
e.g. 48 cm



TuBo-Flex

customised fixation for
endotracheal tubes



Sim-Plex

special size:
e.g. 42 cm



Uni-Fix

mobile DIN rail



Our Team...

We do produce quality – Made in Germany.

Innovation, premium performance and reliability linked with commitment and experience are the statement of our enterprise.

For us customers are not mere final consumers, and our staff are not mere human capital but people who we want to be content and satisfied.

Not only because we learn from our customers' experience thanks to a highly motivated team, but especially because a co-operation coined by mutual trust is more relaxed and just more fun.

„Nothing can sustain someone than the trust offered.“

Paul Claudel

...one Step ahead

*„Dear Ms Schlösser and team,
I would like to express again my sincere thanks for your kindest consultation and the samples you send for free. Great company and really nice people working there.“*
Ms Baston with tracheostomised daughter, 13 month old, Germany

*„Hello,
when I ordered a sample from you about a month ago I was amazed at its prompt and smooth delivery. THANK YOU.
Then I showed your product to my colleagues in the intensive care ward and they were immediately very enthusiastic. Some knew of the product from other hospitals they had previously worked in.
Some were at first rather sceptical that the tubes were fixed firmly enough, but they were astounded at how well the tubes (using the adhesive Velcro® tape) were fixed. It has now been decided to fix tubes using only your holders. The patients now have much fewer or no pressure sores on their mouths and lips. A commendation for your product. THANK YOU“*
Mr. Lemmer, Registered Nurse – Hôpital Kirchberg, Luxembourg

„We tested your samples and found them extremely „nurse-“ and „patient-“ friendly. Beside that they have proven safe and easy to handle. We would definitely like to use your product.“
Dr. Knafelj – University Medical Centre Ljubljana, Slovenia

*„First of all, we used the tracheostomy tube holders with the velcro fastening as the clip type was a little too difficult to attach to the tracheostomy. However we are very pleased with the velcro fasteners and continue to use them today.
The horrendous scars and pressure scores from the previous ties were healed within 10 days of using your product and William has not had another score since.“*
Father of tracheostomised boy William, Great Britain

Always one step ahead – innovative – unique



Hand-Fix, Ute-Fix, BePa-Clip VARIO

Our customers (excerpt from DE)

Klinikum Aschaffenburg
Aschaffenburg, Germany
www.klinikum-aschaffenburg.de

Charité - Universitätsmedizin Berlin
Berlin, Germany
www.charite.de

Universitätsklinikum Bonn
Bonn, Germany
www.ukb.uni-bonn.de

Universitätsklinikum Carl Gustav Carus Dresden
Dresden, Germany
www.uniklinikum-dresden.de

Universitätsmedizin Göttingen Georg-August-Universität
Göttingen, Germany
www.med.uni-goettingen.de

Universitätsklinikum Halle (Saale)
Halle, Germany
www.medizin.uni-halle.de

Medizinische Hochschule Hannover
Hannover, Germany
www.mh-hannover.de

Klinikum der Universität zu Köln
Cologne, Germany
www.uk-koeln.de

Gesundheit Nordhessen Holding AG
Kassel, Germany
www.klinikum-kassel.de

Universitätsmedizin der Johannes Gutenberg-Universität Mainz
Mainz, Germany
www.unimedizin-mainz.de

Klinikum der Universität München Großhadern
Munich, Germany
<http://gch.klinikum.uni-muenchen.de>

Klinikum der Universität Regensburg
Regensburg, Germany
www.uniklinikum-regensburg.de

and others ...

Contact

National Sales

Alte Landstraße 33
50171 Kerpen
Germany
Tel +49 (0)2237 6393-0
Fax +49 (0)2237 6393-30

International Sales / Logistic Center

Hahnenstraße 8
50171 Kerpen
Germany
Tel +49 (0)2237 5623-0
Fax +49 (0)2237 5623-20
E-mail orders@novo-med.de

Headquarter

Alte Landstraße 3
50171 Kerpen
Germany
Tel +49 (0)2237 9287-0
Fax +49 (0)2237 9287-10
E-mail info@novo-med.de

Postal address

Novo Klinik-Service GmbH
P.O. Box 2305
50153 Kerpen
Germany

For sample requests please use our form "Request a sample"
on our homepage
www.novo-med-group.com

