

# Product Focus

April, 2006 Vol. 8

## Keeping you updated...

With our Product Focus Newsletter we hope to keep Biometrix customers exposed to the new products we have to offer, besides our current extensive product line. Our goal is to increase your awareness of our broad range of products in Critical Care, Operating Rooms and Interventional Cardiology.


Medical Field

**Critical  
Care**

## Custom CVC Kits


### Y-Needle

Biometrix Y-Needle Guidewire Introducer unique design allows access to the vein without disconnecting the syringe attached to the needle. The Y-Needle has a one way integral valve to prevent blood leakage and the risk of embolism during insertion.



### CLC 2000® - Positive Pressure Needleless Valve

The **CLC2000®** is a revolutionary connector that provides a constant velocity of fluid which exits the catheter upon disconnection of a luer. This positive displacement feature prevents the backflow of blood into the catheter lumen which can otherwise lead to thrombotic occlusions.



## Contact

### Biometrix B.V.

#### Address

Paardeweide 36  
P.O.B. 7191  
4800 GD Breda  
The Netherlands

#### Phone

31-76-5422720

#### Fax

31-76-5422716

#### E-Mail

[sales@biometrixmedical.com](mailto:sales@biometrixmedical.com)

Visit us on the web! [www.biometrixmedical.com](http://www.biometrixmedical.com)

For more information, samples or a price quotation  
please contact our international sales department at [sales@biometrixmedical.com](mailto:sales@biometrixmedical.com)

4

4 LUMEN CATHETER SETS | ADULT

CODE	DESCRIPTION	Y-NEEDLE	CLC-2000*
HH-8515	8.5Fr x 15cm	-	-
HH-8520	8.5Fr x 20cm	-	-
HH-8530	8.5Fr x 30cm	-	-
HH-I8515	8.5Fr x 15cm	+	-
HH-I8520	8.5Fr x 20cm	+	-
HH-I8530	8.5Fr x 30cm	+	-
HH-E8515	8.5Fr x 15cm	+	+
HH-E8520	8.5Fr x 20cm	+	+
HH-E8530	8.5Fr x 30cm	+	+

3

3 LUMEN CATHETER SETS | ADULT

CODE	DESCRIPTION	Y-NEEDLE	CLC-2000*
HJ-7015	7Fr x 15cm	-	-
HJ-7020	7Fr x 20cm	-	-
HJ-7030	7Fr x 30cm	-	-
HJ-X7015	7.5Fr x 15cm	-	-
HJ-X7020	7.5Fr x 20cm	-	-
HJ-X7030	7.5Fr x 30cm	-	-
HJ-I7015	7.5Fr x 15cm	+	-
HJ-I7020	7.5Fr x 20cm	+	-
HJ-I7030	7.5Fr x 30cm	+	-
HJ-E7015	7.5Fr x 15cm	+	+
HJ-E7020	7.5Fr x 20cm	+	+
HJ-E7030	7.5Fr x 30cm	+	+

3

3 LUMEN CATHETER SETS PEDIATRIC

HJ-5008	5.5Fr x 8cm
HJ-5013	5.5Fr x 13cm
HJ-5030	5.5Fr x 30cm

2

2 LUMEN CATHETER SETS PEDIATRIC

CODE	DESCRIPTION
HK-4005	4Fr x 5cm
HK-4008	4Fr x 8cm
HK-4013	4Fr x 13cm
HK-4030	4Fr x 30cm
HK-5006	5Fr x 6cm
HK-5008	5Fr x 8cm
HK-5013	5Fr x 13cm
HK-5020	5Fr x 20cm
HK-5030	5Fr x 30cm

2

2 LUMEN CATHETER SETS | ADULT

CODE	DESCRIPTION	Y-NEEDLE	CLC-2000*
HK-7015	7Fr x 15cm	-	-
HK-7020	7Fr x 20cm	-	-
HK-7030	7Fr x 30cm	-	-
HK-X7015	7.5Fr x 15cm	-	-
HK-X7020	7.5Fr x 20cm	-	-
HK-X7030	7.5Fr x 30cm	-	-
HK-I7015	7.5Fr x 15cm	+	-
HK-I7020	7.5Fr x 20cm	+	-
HK-I7030	7.5Fr x 30cm	+	-
HK-E7015	7.5Fr x 15cm	+	+
HK-E7020	7.5Fr x 15cm	+	+
HK-E7030	7.5Fr x 30cm	+	+

1

1 LUMEN CATHETER SETS | ADULT

CODE	DESCRIPTION	Y-NEEDLE	CLC-2000*
HO-4015	14Ga x 15cm	-	-
HO-4020	14Ga x 20cm	-	-
HO-4030	14Ga x 30cm	-	-
HO-6015	16Ga x 15cm	-	-
HO-6020	16Ga x 20cm	-	-
HO-6030	16Ga x 30cm	-	-
HO-I4015	14Ga x 15cm	+	-
HO-I4020	14Ga x 20cm	+	-
HO-I4030	14Ga x 30cm	+	-
HO-I6015	16Ga x 15cm	+	-
HO-I6020	16Ga x 20cm	+	-
HO-I6030	16Ga x 30cm	+	-
HO-E4015	14Ga x 15cm	+	+
HO-E4020	14Ga x 20cm	+	+
HO-E4030	14Ga x 30cm	+	+
HO-E6015	16Ga x 15cm	+	+
HO-E6020	16Ga x 20cm	+	+
HO-E6030	16Ga x 30cm	+	+
HO-8015	18Ga x 15cm	-	-
HO-8020	18Ga x 20cm	-	-

1

1 LUMEN CATHETER SETS PEDIATRIC

CODE	DESCRIPTION
HO-9010	20Ga x 10cm
HO-9015	20Ga x 15cm
HO-9020	20Ga x 20cm
HO-7008	24Ga x 8cm
HO-7015	24Ga x 15cm
HO-7020	24Ga x 20cm

1

RAPID INFUSION CATHETER SETS

CODE	DESCRIPTION
JDX-4062	8Fr High flow catheter
JC-7100	Trauma pack

# CENTRA-LINE™


## Central Venous Catheter Sets

## BIOMETRIX

Paardeweide 36 P.O.B. 7191  
4800 GD Breda, The Netherlands  
Tel: 31-76-5422720 | Fax: 31-76-5422716  
sales@biometrixmedical.com  
www.biometrixmedical.com







## Customized CVC Trays

Biometrix offers custom CVC trays including: incise drapes, surgical gowns, gauze balls, gauze swabs, hypodermic needles and needle counters.


# CENTRA-LINE™

## CENTRAL VENOUS CATHETER SETS

Biometrix offers a broad range of catheter kits and trays for central vein insertion. We provide safe and reliable devices for a wide variety of clinical interventions

Wpuncture site:

- CVP monitoring
- Continuous and simultaneous infusion of incompatible drugs
- Blood sampling
- Blood product administration
- TPN administration



## Rapid Infusion Catheter Sets

Rapid Infusion Catheter Sets provide a high flow low rate of fluid and blood. The set is designed for a faster insertion technique intended for rapid central venous catheterization in trauma cases, in emergency rooms or operating theaters.




## CLC 2000®- Positive Pressure Needleless Valve

The CLC 2000® is a revolutionary connector which provides a constant velocity of fluid that exits the catheter upon luer disconnection. This positive displacement feature prevents a back flow of blood into the catheter lumen. Minimizing blood in the lumen, results in less clotting.



## Y-Introducer Needle

The Biometrix Y-Introducer Needle offers a unique design which allows access to the vein without disconnecting the syringe attached to the needle. The Y-Needle has a one way integral valve to prevent the risk of blood leakage and embolism during insertion.



## Kink Resistant Guidewires

Biometrix Kink Resistant Guidewires offer considerable advantages. Designed to withstand sharp angles without buckling and coiling, thus facilitating passage of the CVC through the patient's vasculature.