

# Cover-ENT<sup>®</sup>

Disposable sterile sheaths for flexible nasopharyngoscopes

Sterile, practical and ready to use


Enhances safety & productivity

Preserves efficiency of the nasopharyngoscopes

Reduces service expenses and losses of time

Effective barrier against viruses and bacteriums

Less risk of infection and adverse events



# Cover-ENT® Disposable sterile sheaths for flexible nasopharyngoscopes

*We can provide custom specific models upon customer's requests*

Flexible Endoscopes	Endoscope code	Cover sheath code
<b>OLYMPUS</b>	ENF-XP	GNF20X30
	ENF-P4	GNF42X30O
	ENF-GP	GNF42X30O
	ENF-GP2	GNF33X30O
	ENF-V	GNF44X37
	ENF-V2	GNF42X30O
	ENF-V3	GNF30X30P
	ENF-V4	GNF30X32O
	ENF-VH	GNF4230L2
	ENF-VQ	GNF4230L3 / GNF4430L2
<b>STORZ</b>	11101SK-SK2	GNF30X27
	11101RP1-2	GNF39X30S
	11101CM	GNF39X30
	11101HD	GNF39X30
	11101VP	GNF39X30V
	11102CM	GNF30X30

Flexible Endoscopes	Endoscope code	Cover sheath code	
<b>PENTAX</b>	FNL-7RP3	GNF30X30S	
	FNL-10RP3	GNF42X30O	
	FNL-10RBS	GNF42X30O	
	VNL1170K	GNF4430L2	
	VNL1190K	GNF4430L2	
<b>XION</b>	EF-N SLIM	GNF30X32	
	EF-N / EF-NS	GNF39X32	
	EV-N	GNF44X32	
	EV-NC	GNL39X30S	
	<b>FENTEX</b>	700710FX	GNF35X30
		700761FX	GNF30X30P
700721FX Video		GNF39X30	
<b>MACHIDA</b>	ENT-30PIII	GNF35X30M	
<b>OTOPRONT</b>	VN-S	GNF39X28	
<b>ORL VISION</b>	RS1-RX1	GNFORLV01	
<b>OPTIM</b>	Entity	GNF39X30E	
<b>HENKE SASS WOLF</b>	8300019728	GNF35X30	

## Packaging:

*A dispenser contains 20 individually packaged sheaths, each one with double stickers for test traceability. A box contains 5 dispensers.*

Produced by:

**COVER**

Via Nizza, 102 - 10126 Torino (Italy)  
Tel +39 0182 077213 - Fax +39 0182 989819  
email: info@cover-srl.com web site: www.cover-srl.com

**SINCERT**



**CE**  
0476

UNI CEI EN ISO  
13485:2012

Distributed by: