

Preventing healthcare
associated infections
with our filter portfolio

Our new generation of breathing system filters

Quality at every corner

QUALITY & SAFETY



- **Quality made in Germany**
- **Manufactured at Dräger in Lübeck**
 - 100 % test for leakages
 - 100 % tests of retention rates
 - 100 % test of resistance
- **No spot checks, individual control of every single filter**

SUSTAINABILITY



- **CO₂-footprint improved by:**
 - reduced weight
 - short delivery routes
 - optimised packaging
 - waste-optimised design





PATIENT SAFETY



- **Reduced weight, thus:**
 - low risk of tube dislocation
- **Transparent housing for better visual inspection**

PRODUCT ADVANTAGES



- **Writable surface to document time of filter application**
- **Highly effective bacterial and viral retention rates**
- **Clear visibility of filter type, dead-space, single patient use**
- **Lower moisture loss with HME**



Here you will find the new part numbers

Previous part number	Previous produkt name	New product name	New part number
MP01730	HME HumidStar 55	HME HumidStar® 55 Plus, ~55 ml deadspace	MP05730
MP01735	HME HumidStar 25	HME HumidStar® 25 Plus, ~25 ml deadspace	MP05735
MP01755	Filter CareStar 45	Electrostatic Filter CareStar® 35 Plus, ~35 ml deadspace	MP05755
MP01770	Filter CareStar 30	Electrostatic Filter CareStar® 20 Plus, ~20 ml deadspace	MP05770
MP01785	Filter SafeStar 80	Mechanical Filter SafeStar® 90 Plus, ~90 ml deadspace	MP05785
MP01790	Filter SafeStar 55	Mechanical Filter SafeStar® 55 Plus, ~55 ml deadspace	MP05790
MP01795	Filter SafeStar 60A	Mechanical Filter Safe Star® 60A Plus, ~60 ml deadspace, angled	MP05795
MP01800	Filter/HME TwinStar 90	Electrostatic Filter and HME TwinStar® 90 Plus, ~90 ml deadspace	MP05800
MP01801	Filter/HME TwinStar HEPA	Mechanical Filter and HME TwinStar® HEPA Plus, ~90 ml deadspace	MP05801
MP01805	Filter/HME TwinStar 55	Electrostatic Filter and HME TwinStar® 55 Plus, ~55 ml deadspace	MP05805
MP01810	Filter/HME TwinStar 65A	Electrostatic Filter and HME TwinStar® 60A Plus, ~60 ml deadspace, angled	MP05810
MP01815	Filter/HME TwinStar 25	Electrostatic Filter and HME TwinStar® 25 Plus, ~25 ml deadspace	MP05815
MP01820	Filter/HME TwinStar 8	Electrostatic Filter and HME TwinStar® 9 Plus, ~9 ml deadspace	MP05820

Packaging unit = 100 Pcs.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Locate your Regional Sales Representative at:
www.draeger.com/contact



REGION EUROPE

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel. +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel. +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA PACIFIC

Dräger Singapore Pte. Ltd.
61 Science Park Road
The Galen #04-01
Singapore 117525
Tel. +65 6872 9288
Fax +65 6259 0398
asia.pacific@draeger.com

REGION CENTRAL

AND SOUTH AMERICA
Dräger Indústria e Comércio Ltda.
Al. Pucuruí - 51 - Tamboré
06460-100 - Barueri - Sao Paulo
Tel. +55 (11) 4689-4900
relacionamento@draeger.com