

PROFEEL[®]

DHD[™] Micro Powder Free
Latex Surgical Gloves

Powder Free Latex Sterile Surgical Gloves

PROFEEL Powder Free

PROFEEL DHD[™] Platinum

PROFEEL DHD[™] Micro

PROFEEL DHD[™] Extra Protection

PROFEEL DHD[™] Sensitive



Damp Hand Donning – DHD[™]

Natural Rubber Latex

Powder Free



Your Partner In Protection[™]

Extra Sensitivity

The approach taken in meeting the healthcare professionals' increasingly stringent quality requirement is reflected in the outstanding features of **PROFEEL® DHD™ Micro Powder Free Latex Surgical Gloves**.

Thinner than normal surgical gloves, **PROFEEL® DHD™ Micro Powder Free Latex Surgical Gloves** are specially designed for specialists who demand extra sensitivity in the fingertips for micro and ophthalmic surgeries. The optimum strength-to-thickness ratio ensures that the basic barrier protection property is not compromised. The unique brown colour minimizes eye-strain in the Operating Theatre, provides a good contrast to light-coloured organs and tissues and aids visibility in surgical works

- Non-slick, skin grip surface finish provides good instrument grip

- Thinner than normal surgical gloves for added sensitivity

- High tensile strength minimizes tear

- Anatomically shaped to reduce hand fatigue during prolonged surgical work

- Printed size and logo for easy identification

- Reinforced cuff prevents tearing while donning

- Beaded cuff prevents tearing while donning

PRODUCT INFORMATION

Type: Powder free, curved fingers, bisque finish and beaded cuff

Primary Material: Natural rubber latex

Dimensions:

Glove Size	Reorder #	Palm Width (mm)	Length (mm)
5½	P3355-34	72 ± 4	278 min
6	P3360-34	77 ± 5	280 min
6½	P3365-34	83 ± 5	280 min
7	P3370-34	89 ± 5	283 min
7½	P3375-34	95 ± 5	287 min
8	P3380-34	102 ± 6	288 min
8½	P3385-34	108 ± 6	290 min
9	P3390-34	114 ± 6	290 min

Thickness:

Location of Thickness Measurements	Single Wall (mm)
Finger (at 15 mm from the extreme tip)	0.18 ± 0.03
Palm (at center of palm)	0.15 ± 0.03
Cuff (at 25 mm from cuff end)	0.13 ± 0.03

Typical Physical Properties:

Parameters	Before Aging	After Aging
Tensile Strength (MPa)	29	27
Force At Break (N)	15	14

Powder: Polymer coating online and offline (DHD-Micro-Brown)

Color: Brown (PMS 4635)

Design and Feature: Hand specific, curved fingers, bisque finish and beaded cuff

Packing: 1 pair (1 left and 1 right) of gloves per inner wrap, 1 inner wrap per pouch, 50 pouches per dispenser, 4 dispensers per carton

Pre-shipment Quality Inspection:

Criteria	Sampling Plan
Dimension	N=13 (EN455-2); S-2, AQL 4.0
Physical Properties	N=13 (EN455-2); S-2, AQL 4.0
1000ml Water Leak Test	G-I, AQL 1.0
Protein Content	N=8 (EN455-3)
Powder Free Residue	N=3 pairs
Major Visual Inspection	G-I, AQL 1.5
Minor Visual Inspection	G-I, AQL 4.0

Product Conformance : • Personal Protective Equipment of Complex Design Category III, in compliance with 89/686/EEC, type tested to EN 420:2003+A1:2009, EN 374-2:2003, EN 374-3:2003 & EN 388:2003, CE0086

• Medical Device, in compliance with 93/42/EEC, type tested to EN 455-1:2000, EN 455-2:2009, EN 455-3:2006, EN 455-4:2009 and ASTM D3577, CE0086 Class IIa Sterile

Quality Assurance : • US FDA Quality System Regulation (QSR)
• ISO9001 Quality Management System
• ISO13485 Quality Management System



BS EN ISO 9001:2000
BRITISH STANDARDS INSTITUTION
CERTIFICATE NO. FM 13934



REGISTERED TO BRITISH STANDARD
BS EN ISO 13485:2012 PARTS 1, 2 & 3



MKS BPP, UK



GKL
BS EN 374-3:2003
Chemical Hazards



BS EN 374-2:2003
Micro-organism
Hazards



BS EN 388:2003
Mechanical
Hazards



WRP Asia Pacific Sdn Bhd 1 4 7 8 1 7 V

Lot 1, Jalan 3, Kawasan Perusahaan
Bandar Baru Salak Tinggi,
43900 Sepang, Selangor Darul Ehsan,
MALAYSIA

Office +60-3-8706 1486

Facsimile +60-3-8706 1557

Email customer_wrp@wrpworld.com



REMESCO Handelges.m.b.H.

Grinzinger Allee 5/25, 1190 Vienna,
AUSTRIA

Office +43-1-328-50-88

Facsimile +43-1-328 29-49

Email office@remesco.com

Distributed by: