



MEDIBORGH



Reusable isolatiejas

Item : DR50009-MS

Specifications

- Reusable, 75-100 times washable
- AAMI PB70 Rated Level 1 / Level 2
- DWR Finishing, >Point70 (Grade 3)
- 99% polyester / 1% carbon with inherent static control properties / 5mm stripes Overlap back, twill tie neck and back waist closures
- White polyester knit cuffs

A high-count, densely woven fabric made with non-linting continuous filament polyester fibers that offer excellent fluid resistance.

This fabric is also engineered with carbon fibers to offer protection where static control is required.



Instructions for washing and sterilisation

Wash with warm water, the maximum temperature shall not exceed 80 °C

During **disinfection**, it is necessary to avoid strong alkaline **disinfectants**, which will damage the fabric. It is recommended to use ethylene oxide sterilization, steam sterilization or irradiation sterilization.



Binnendienst@Mediborgh.nl



0031 74 349 25 22



www.mediborgh.nl