

Sharpak COMPACT, COMPACT 'O', 6, 18 PLUS and 18 ECO sharps containers

The Sharpak range of blue-top sharps containers are recommended for the segregation, collection and disposal of partially discharged sharps waste which are NOT contaminated with cytotoxic or cytostatic medicines.

This comprehensive range of 10 containers is available from 0.5L through to 26.3L and is accompanied by a range of brackets, trays and stands.

Sharpak sharps containers are approved to ISO 23907:2012 and tested and certified in accordance with the provisions of the United Nations recommendations on the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified, NOS.

Sharpak sharps containers MUST only be disposed of by incineration.

Please note:- For disposal of sharps contaminated with cytotoxic or cytostatic medicines please refer to the data sheet for the Code PURPLE range of Sharpak HYDRI sharps containers.



Patented and trademarked by Nolato.

MANUFACTURED, TESTED AND CERTIFIED TO ISO 23907:2012

UN TYPE TESTED

BLUE COLOUR CODED CLOSURES TO MEET THE RECOMMENDATIONS ON WASTE SEGREGATION AS DETAILED IN THE DEPARTMENT OF HEALTH'S TECHNICAL MEMORANDUM HTM-07-01

SECURE CLOSURE WITH LOCKING FACILITY

VISIOLOCK INDICATOR IS FITTED ON THE SHARPAK COMPACT, COMPACT 'O', 18 PLUS AND 18 ECO

DISTINCTIVE COLOURING AND CLEAR LABELLING

SEQUENTIAL NUMBERING FOR TRACEABILITY

TOUGH AND DURABLE

SUPPLIED ASSEMBLED

COST EFFECTIVE

PURPOSE DESIGNED WALL/TROLLEY BRACKETS, COMPACT STAND, PROCEDURE TRAYS AND WALL UNIT

INSULIN PEN NEEDLE REMOVING FEATURE ON SHARPAK 18 ECO

ALSO AVAILABLE WITHOUT NON RETURN PETALS (OPEN APERTURE) COMPACT 'O'



0.5 LITRES

0.6 LITRES





1.53 LITRES

1.58 LITRES



PRODUCT	MAX CAPACITY	CAPACITY TO FILL LINE	PRODUCT CODE	QTY	HEIGHT	WIDTH	DEPTH
Sharpak COMPACT	0.5 litres	0.39 litres	SU06201	24	186 mm	67 mm	67 mm
Sharpak COMPACT 'O'	0.5 litres	0.39 litres	SU06202	24	186 mm	67 mm	67 mm
Sharpak 6	0.6 litres	0.48 litres	SU07400	60	200 mm	102 mm	46 mm
Sharpak 18 ECO	1.53 litres	1.22 litres	SU06496	24	190 mm	104 mm	104 mm
Sharpak 18 PLUS	1.58 litres	1.27 litres	SU06495	48	210 mm	104 mm	104 mm

ACCESSORIES	CODE	QTY	ACCESSORIES	CODE	QTY
Wall & trolley bracket for 18 PLUS	SM06500	1	Compact Stand	SM07500	1
Shaped tray for 18 PLUS & COMPACT	SU06211	5	Sharps Injury Prevention Pack (SIPP)	SB06196	6
SAFE & SOUND wall unit 18 PLUS	SS00018	1			



Actual aperture sizes when printed at 100%

Dimensions in mm

For details regarding Sharpak distributors and worldwide availability, please contact:
Nolato Jaycare Ltd Walton Road Portsmouth PO6 1TS
Tel +44 2392 314711 Fax +44 2392 380314 www.sharpaks.com