


# Structure >>



## Transition joint >>



## Product Introduction >>

- Unique auto-inflation device enables the operator to inflate the cuff with one hand.
- Meet the needs of lung isolation in patients of different ages.
- Sealed fiber bronchoscope port for insertion into a fiber bronchoscope or suction tube.
- Reliable endotracheal tube mount to avoid occlusion displacement during surgery.
- Flexible and slim design of the main tube to adjust the position and orientation of the blocker according to your needs.
- Innovative Tappa angel design, in line with physiological anatomy, easier for intubation.
- High volume low pressure silicone cuff to reduce bronchial mucosal irritation and damage.
- Used with endotracheal tube to avoid secondary damage caused by tube replacement.

